



Company Name: Maymom LLC
Postal Address: 20-7, Aly 19, LN 333, Chungcheng Rd
Chunan Township, Miaoli County Taiwan
Zip Code: 350
Telephone number: +886 37470556
Email: CS@maymom.com

Cleaning Statement

Date of issue: Apr 23, 2025

To whom it may concern:

We, the manufacturer, declare that our flange products and the items made of polypropylene and silicone mentioned below can be sterilized in the microwave steam bags, boiling water or autoclave 121 degrees Celsius for 15 minutes.

Product SKU (Model Number):

Polypropylene (PP) Products:

A012-Fxx, A012-B, L012W-B, MP012C-Fxx, MP012-Fxx, M012-B, M012-tB, M012-Fxx, S012W-Fxx, S012W-B, S012W.2G-B

Note: "xx" represents the size of the flanges. Packs of multiple units are not listed, but their nomenclature follows the same prefix.

Silicone Products:

M003-M, S005.40-D, A005-D, E005-D, L005-D, P005.2G-D, XP024-M

Polypropylene (PP) and Silicone Products:

M005-2V4M, M005-4V6M

We have done autoclave study on our flanges/shields as the report below.

Reference: report #. TWNC00596197

URL: <http://www.maymom.com/en/reports/TWNC00596197%20autoclave.pdf>

All the polypropylene products mentioned above use the same material so they can fall into the same category as what is on the report.

FDA-compliant reports for our silicone-based products indicate a significantly higher melting point compared to polypropylene. Since the product launch, we have not received any customer feedback regarding damage to the silicone products from steam-based sterilization. Therefore, we recommend using the same autoclave conditions for these products.

Be sure NOT to set the autoclave temperature/pressure above what is indicated in the report. Use 121 degrees Celsius, 14.2 psi, 15 min for the sterilization process.

Name: Liangchi Li

Function: Manger of Maymom LLC

Signature: